

Collagen III alpha 1/COL3A1 Rabbit mAb

Catalog No.: A0817 **Recombinant** **33 Publications**

Basic Information

Observed MW

150 kDa

Calculated MW

139 kDa

Category

Primary antibody

Applications

WB,Auto WB,IP,ELISA

Cross-Reactivity

Human, Mouse, Rat

CloneNo number

ARC2542

Background

This gene encodes the pro-alpha1 chains of type III collagen, a fibrillar collagen that is found in extensible connective tissues such as skin, lung, uterus, intestine and the vascular system, frequently in association with type I collagen. Mutations in this gene are associated with Ehlers-Danlos syndrome type IV, and with aortic and arterial aneurysms.

Recommended Dilutions

WB 1:1000 - 1:6000

Auto WB 1:100 - 1:500

IP 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells

ELISA Recommended starting
concentration is 1 µg/mL.
Please optimize the
concentration based on
your specific assay
requirements.

Immunogen Information

Gene ID

1281

Swiss Prot

P02461

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

EDS4A; EDSVASC; PMGEDSV; Collagen III alpha 1/COL3A1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide, 0.05% BSA, 50% glycerol, pH7.3.


Validation Data



Immunoprecipitation analysis of 300 µg extracts from HepG2 cells using 3 µg Collagen III alpha 1/COL3A1 Rabbit mAb (A0817). Western blot was performed from the immunoprecipitate using Collagen III alpha 1/COL3A1 Rabbit mAb (A0817) at a dilution of 1:1000.




Western blot analysis of various lysates using Collagen III alpha 1/COL3A1 Rabbit mAb (A0817) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s.



Western blot analysis of various lysates using Collagen III alpha 1/COL3A1 Rabbit mAb (A0817) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 60s.

Validation Data



Simple Western™ analysis of lysates from Hep G2 cells using Collagen III alpha 1/COL3A1 Rabbit mAb (A0817) at 1:100 dilution. The virtual lane view (left) shows the target band (as indicated) with samples in concentrations of 0.2 mg/mL, 0.7 mg/mL and 2.0 mg/mL. The corresponding electropherogram view (right) plots chemiluminescence intensity against molecular weight along the capillary for sample concentrations of 0.2 mg/mL, 0.7 mg/mL and 2.0 mg/mL. This experiment was performed under reducing conditions on the Jess™ Simple Western instrument from ProteinSimple, a BioTechne brand, using the 12-230 kDa separation module.