

GNG12 Rabbit pAb

Catalog No.: A10651 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

8kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables PDZ domain binding activity. Predicted to be involved in G protein-coupled receptor signaling pathway. Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

55970

Swiss Prot

Q9UBI6

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HG3A; GNG12

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.