

AEN Rabbit pAb

Catalog No.: A15182

Basic Information

Observed MW

Refer to figures

Calculated MW

36kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables exonuclease activity. Involved in intrinsic apoptotic signaling pathway in response to DNA damage by p53 class mediator and response to ionizing radiation. Located in nuclear membrane; nucleolus; and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

64782

Swiss Prot

Q8WTP8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ISG20L1; pp12744; AEN

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.