

SHC3 Rabbit pAb

Catalog No.: A15453

Basic Information

Observed MW

Refer to figures

Calculated MW

64kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

Enables phosphotyrosine residue binding activity. Predicted to be involved in transmembrane receptor protein tyrosine kinase signaling pathway. Predicted to act upstream of or within glutamatergic synaptic transmission and learning or memory. Predicted to be located in cytosol. Predicted to be active in plasma membrane.

Recommended Dilutions

WB 1:200 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

53358

Swiss Prot

Q92529

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

RAI; NSHC; SHCC; N-Shc; SHC3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.