

PE Mouse anti-Human CD269(BCMA) mAb

Catalog No.: A22844

Basic Information

Observed MW

Refer to figures

Calculated MW

20kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

AMC0550

Conjugate

PE. Ex:565nm. Em:574nm.

Recommended Dilutions

FC 10 µL/Test

Background

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is preferentially expressed in mature B lymphocytes, and may be important for B cell development and autoimmune response. This receptor has been shown to specifically bind to the tumor necrosis factor (ligand) superfamily, member 13b (TNFSF13B/TALL-1/BAFF), and to lead to NF-kappaB and MAPK8/JNK activation. This receptor also binds to various TRAF family members, and thus may transduce signals for cell survival and proliferation.

Immunogen Information

Gene ID

608

Swiss Prot

Q02223

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

BCM; BCMA; CD269; TNFRSF13A

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

IgG1

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.1% Proclin300, 0.2% BSA, pH7.2.

Validation Data



Flow cytometry: 1X10⁶ Human PBMC were surface-stained with PE Mouse anti-Human CD269(BCMA) mAb (A22844, orange line) or PE Mouse IgG1 isotype control mAb (A22878, blue line).