

PE Rabbit anti-Mouse CD130/gp130 mAb

Catalog No.: A25264 **1 Publications**

Basic Information

Observed MW

Calculated MW
102kDa**Category**
Primary antibody**Applications**
FC**Cross-Reactivity**
Mouse**CloneNo number**
ARC64569**Conjugate**
PE. Ex:565nm. Em:574nm.

Background

Contributes to interleukin-11 binding activity and interleukin-11 receptor activity. Involved in interleukin-6-mediated signaling pathway; intestinal epithelial cell development; and positive regulation of Notch signaling pathway. Acts upstream of or within several processes, including glycogen metabolic process; positive regulation of astrocyte differentiation; and regulation of Notch signaling pathway. Located in dendrite; external side of plasma membrane; and neuronal cell body. Is expressed in several structures, including alimentary system; brain; early conceptus; immune system; and reproductive system. Used to study rheumatoid arthritis and stomach cancer. Human ortholog(s) of this gene implicated in arteriosclerosis; hyper IgE recurrent infection syndrome 4; and hyperandrogenism. Orthologous to human IL6ST (interleukin 6 cytokine family signal transducer).

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Immunogen Information

Gene ID
16195**Swiss Prot**
Q00560

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CD130; gp130; D13Ertd699e; 5133400A03Rik

Contact

 | 400-999-6126
 | cn.market@abclonal.com.cn
 | www.abclonal.com.cn

Product Information


Source
Rabbit**Isotype**
IgG**Purification**
Affinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Flow cytometry: 1×10^6 Neuro-2a cells (Low Expression, left) and NIH/3T3 cells (right) were surface-stained with PE Rabbit anti-Mouse CD130/gp130 mAb (A25264, 5 μ l/Test, orange line) or PE Rabbit IgG isotype control (A24172, 5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1×10^6 NIH/3T3 cells were surface-stained with PE Rabbit IgG isotype control (A24172, 5 μ l/Test, left) or PE Rabbit anti-Mouse CD130/gp130 mAb (A25264, 5 μ l/Test, right).