

XPB1 Rabbit mAb

Catalog No.: A25319

Recombinant

2 Publications

Basic Information

Observed MW

60kDa/56kDa/29kDa,60kDa

Calculated MW

29kDa

Category

Primary antibody

Applications

WB,IP,ELISA

Cross-Reactivity

Human, Mouse

CloneNo number

ARC3254

Background

This gene encodes a transcription factor that regulates MHC class II genes by binding to a promoter element referred to as an X box. This gene product is a bZIP protein, which was also identified as a cellular transcription factor that binds to an enhancer in the promoter of the T cell leukemia virus type 1 promoter. It may increase expression of viral proteins by acting as the DNA binding partner of a viral transactivator. It has been found that upon accumulation of unfolded proteins in the endoplasmic reticulum (ER), the mRNA of this gene is processed to an active form by an unconventional splicing mechanism that is mediated by the endonuclease inositol-requiring enzyme 1 (IRE1). The resulting loss of 26 nt from the spliced mRNA causes a frame-shift and an isoform XPB1(S), which is the functionally active transcription factor. The isoform encoded by the unspliced mRNA, XPB1(U), is constitutively expressed, and thought to function as a negative feedback regulator of XPB1(S), which shuts off transcription of target genes during the recovery phase of ER stress. A pseudogene of XPB1 has been identified and localized to chromosome 5.

Recommended Dilutions

WB 1:500 - 1:1000**IP** 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells**ELISA** Recommended starting
concentration is 1 µg/mL.
Please optimize the
concentration based on
your specific assay
requirements.

Immunogen Information

Gene ID

7494

Swiss Prot

P17861

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

XPB2; TREB5; XPB-1; TREB-5; XPB1s

Contact

 | 400-999-6126 | cn.market@abclonal.com.cn | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.


Validation Data



Immunoprecipitation of XBP1 from 200 µg extracts of Huh7 cells treated with thapsigargin (300 nM, 18h) was performed using 0.5 µg of XBP1 Rabbit mAb (A25319). Rabbit IgG isotype control (AC042) was used to precipitate the Control IgG sample. IP samples were eluted with 1X non-reducing Laemmli Buffer. The Input lane represents 10% of the total input. Western blot analysis of immunoprecipitates was conducted using XBP1 Rabbit mAb (A25319) at a dilution of 1:1000.



Western blot analysis of lysates from HeLa cells using XBP1 Rabbit mAb (A25319) at 1:1000 dilution incubated overnight at 4°C. HeLa cells were treated with tunicamycin (10 µg/ml) for 8 hours. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 30 µg per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 60 s.



Western blot analysis of lysates from RAW 264.7 cells using XBP1 Rabbit mAb (A25319) at 1:1000 dilution incubated overnight at 4°C. RAW264.7 cells were treated with tunicamycin (5 µg/ml) for 8 hours. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 30 µg per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 90s.