

ABflo® 594 Mouse IgG 1 isotype control

Catalog No.: A25489

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

CloneNo number

AMC50005

Conjugate

ABflo® 594. Ex:588nm. Em:604nm.

Background

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300, pH7.3.

Validation Data



Flow cytometry: 1×10^6 MCF7 cells (negative control, left) and SH-SY5Y cells (right) were surface-stained with ABflo® 594 Mouse Anti-Human CD56 mAb (5 µl/Test, orange line) or ABflo® 594 Mouse IgG 1 isotype control (A25489, 5 µl/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1×10^6 SH-SY5Y cells were surface-stained with ABflo® 594 Mouse IgG 1 isotype control (A25489, 5 µl/Test, left) or ABflo® 594 Mouse Anti-Human CD56 mAb (5 µl/Test, right).