

ZNF217 Rabbit pAb

Catalog No.: A27151

Basic Information

Observed MW

Calculated MW
115kDa

Category
Primary antibody

Applications
WB, ELISA

Cross-Reactivity
Human

Background

Enables DNA-binding transcription repressor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Involved in negative regulation of transcription by RNA polymerase II. Located in mitochondrion and nuclear speck. Part of histone deacetylase complex.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID	Swiss Prot
7764	075362

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms
ZABC1

Contact

	400-999-6126
	cn.market@abclonal.com.cn
	www.abclonal.com.cn

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.09% Sodium azide, 50% glycerol, pH 7.3.