

APC Rabbit anti-Human IL-17A mAb

Catalog No.: A27825

Basic Information

Observed MW**Calculated MW**

18kDa

Category

Primary antibody

Applications

FC (intra)

Cross-Reactivity

Human

CloneNo number

ARC70526

Conjugate

APC. Ex:650nm. Em:660nm.

Background

This gene is a member of the IL-17 receptor family which includes five members (IL-17RA-E) and the encoded protein is a proinflammatory cytokine produced by activated T cells. IL-17A-mediated downstream pathways induce the production of inflammatory molecules, chemokines, antimicrobial peptides, and remodeling proteins. The encoded protein elicits crucial impacts on host defense, cell trafficking, immune modulation, and tissue repair, with a key role in the induction of innate immune defenses. This cytokine stimulates non-hematopoietic cells and promotes chemokine production thereby attracting myeloid cells to inflammatory sites. This cytokine also regulates the activities of NF-kappaB and mitogen-activated protein kinases and can stimulate the expression of IL6 and cyclooxygenase-2 (PTGS2/COX-2), as well as enhance the production of nitric oxide (NO). IL-17A plays a pivotal role in various infectious diseases, inflammatory and autoimmune disorders, and cancer. High levels of this cytokine are associated with several chronic inflammatory diseases including rheumatoid arthritis, psoriasis and multiple sclerosis. The lung damage induced by the severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is to a large extent, a result of the inflammatory response promoted by cytokines such as IL17A.

Recommended Dilutions

FC (intra) 5 µl per 10⁶ cells in
100 µl volume

Immunogen Information

Gene ID

3605

Swiss Prot

Q16552

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

IL17; CTLA8; IL-17; ILA17; CTLA-8; IL-17A

Contact

 | 400-999-6126 | cn.market@abclonal.com.cn | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data



Flow cytometry: 1X10⁶ Human PBMC (untreated,left) and Human PBMC (50ng/mL PMA and 1ug/mL ionomycin for 4-6 hours in the presence of monensin,right) were surface-stained with APC Rabbit anti-Human IL-17A mAb (A27825,5 µl/Test)

Flow cytometry: 1X10⁶ Human PBMC (untreated,left) and Human PBMC (treated with 50 ng/ml PMA and 1 µg/ml ionomycin for 6 hours and 2µM Monensin for last 5 hour,right) were intracellularly-stained with ABflo® 488 Rabbit anti-Human IL-2 mAb (A27825,0.5 µg).