

PE Rabbit anti-NHP IgD mAb

Catalog No.: A28180

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

Cynomolgus, Rhesus

CloneNo number

ARC60233

Conjugate

PE. Ex:565nm. Em:574nm.

Background

Recommended Dilutions

FC ≤0.25 µg per million cells
 in 100 µl volume

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.


Validation Data



Flow cytometry: 1×10^6 Cynomolgus PBMC were surface-stained with APC Rabbit anti-Human/Monkey CD19 mAb (A26220, 5 μ l/Test) and PE Rabbit IgG isotype control (A24172, 5 μ l/Test, left) or PE Rabbit anti-NHP IgD mAb (A28180, 0.25 μ g, right). Cells in the lymphocyte gate were used for analysis.



Flow cytometry: 1×10^6 Rhesus PBMC were surface-stained with ABflo® 450 Mouse anti-Human CD11b mAb (A27475, 5 μ l/Test), ABflo® 488 Rabbit anti-Human/Monkey CD19 mAb (A23008, 5 μ l/Test) and PE Rabbit IgG isotype control (A24172, 5 μ l/Test, left) or PE Rabbit anti-NHP IgD mAb (A28180, 0.25 μ g, right). Cells in the CD11b- gate were used for analysis.



Flow cytometry: 1×10^6 Human peripheral blood mononuclear cells (negative control, left) and Cynomolgus peripheral blood mononuclear cells (right) were surface-stained with PE Rabbit anti-NHP IgD mAb (A28180, 0.25 μ g, right). Cells in the lymphocyte gate were used for analysis.