

PE Rabbit anti-Mouse CD117/c-Kit mAb

Catalog No.: A28883

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ARC60551

Conjugate

PE. Ex:565nm. Em:574nm.

Background

The c-Kit proto-oncogene is the cellular homolog of the transforming gene of a feline retrovirus (v-Kit). The c-kit protein includes characteristics of a protein kinase transmembrane receptor. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

FC ≤0.25 µg per million cells
in 100 µl volume

Immunogen Information

Gene ID

16590

Swiss Prot

P05532

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

W; Bs; Fdc; Ssm; SC01; SC05; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfsc03

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data



Flow cytometry: 1X10⁶ C57BL/6 mouse bone marrow cells were surface-stained with ABflo® 450 Rat anti-Mouse CD45R/B220 mAb (A27783, 0.5 µg) and PE Rabbit IgG isotype control (A24172, 5 µl/Test, left) or PE Rabbit anti-Mouse CD117/c-Kit mAb (A28883, 0.25 µg, right). Total viable cells were used for analysis.