

Zinc- α 2-glycoprotein (ZAG/AZGP1) Rabbit pAb

Catalog No.: A5365

1 Publications

Basic Information

Observed MW

41kDa

Calculated MW

34kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Involved in cell adhesion and detection of chemical stimulus involved in sensory perception of bitter taste. Located in extracellular space.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

563

Swiss Prot

P25311

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ZAG; ZA2G; Zinc- α 2-glycoprotein (ZAG/AZGP1)

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of lysates from human serum, using Zinc-α2-glycoprotein (ZAG/AZGP1) Rabbit pAb (A5365) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s.