

MAPK15 Rabbit pAb

Catalog No.: A8286

Basic Information

Observed MW

Refer to figures

Calculated MW

60kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables MAP kinase activity and chromatin binding activity. Involved in several processes, including dopamine uptake; positive regulation of DNA metabolic process; and regulation of organelle organization. Located in several cellular components, including Golgi apparatus; autophagosome; and microtubule organizing center. Biomarker of breast cancer.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

225689

Swiss Prot

Q8TD08

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ERK7; ERK8; MAPK15

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.