

# Mouse anti HA-Tag mAb

Catalog No.: AE008

269 Publications

## Basic Information

### Observed MW

26 kDa

### Calculated MW

### Category

Tag antibody

### Applications

WB,IP,IF/ICC,ELISA

### Cross-Reactivity

Species independent

### CloneNo number

AMC0503

## Background

Protein tags are peptide sequences genetically grafted onto a recombinant protein. Often these tags are removable by chemical agents or by enzymatic means, such as proteolysis or intein splicing. Tags are attached to proteins for various purposes. Epitope tags are short peptide sequences which are chosen because high-affinity antibodies can be reliably produced in many different species. These are usually derived from viral genes, which explain their high immunoreactivity. Epitope tags include V5-tag, Myc-tag, HA-tag and NE-tag. These tags are particularly useful for western blotting, immunofluorescence and immunoprecipitation experiments, although they also find use in antibody purification.

## Recommended Dilutions

**WB** 1:10000 - 1:20000

**IP** 0.5 µg-4 µg antibody for  
200 µg-400 µg extracts  
of whole cells

**IF/ICC** 1:50 - 1:200

**ELISA** Recommended starting  
concentration is 1 µg/mL.  
Please optimize the  
concentration based on  
your specific assay  
requirements.

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Immunogen Information

### Gene ID

### Swiss Prot

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

HA;HA tag;HA-tag

## Product Information

### Source

Mouse

### Isotype

IgG1,Kappa

### Purification


Affinity purification

### Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.


## Validation Data




Western blot analysis of over-expressed HA-tagged protein in 293T cell using Mouse anti HA-Tag mAb (AE008) at different dilution. Each lane was loaded with 2 ug cell lysate.  
 Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS003) at 1:10000 dilution.  
 Blocking buffer: 3% nonfat dry milk in TBST.  
 Detection: ECL Basic Kit (RM00020).  
 Exposure time: 1s.



Western blot analysis of lysates from 293T cells, using Mouse anti HA-Tag mAb (AE008) at 1:10000 dilution.  
 Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS003) at 1:10000 dilution.  
 Lysates/proteins: 25µg per lane.  
 Blocking buffer: 3% nonfat dry milk in TBST.  
 Detection: ECL Basic Kit (RM00020).  
 Exposure time: 10s.



Immunofluorescence analysis of 293T cells transfected with HA-Tag fusion protein and untreated 293T cells use Mouse anti HA-Tag mAb (AE008) at dilution of 1:50 (40x lens).  
 Blue: DAPI for nuclear staining.



Immunoprecipitation of over-expressed HA-tagged protein in 293T cells using HA-tag antibody (AE008). A mock served as negative control and over-expressed 293T cell lysate served as positive control.