

Phospho-IRS1-S307 Rabbit pAb

Catalog No.: AP0371 **2 Publications**

Basic Information

Observed MW

Calculated MW
132kDa

Category
Primary antibody

Applications
ELISA

Cross-Reactivity
Human, Mouse, Rat

Background

This gene encodes a protein which is phosphorylated by insulin receptor tyrosine kinase. Mutations in this gene are associated with type II diabetes and susceptibility to insulin resistance.

Recommended Dilutions

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID 3667	Swiss Prot P35568
------------------------	-----------------------------

Immunogen
Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms
HIRS-1; Phospho-IRS1-S307

Contact

		400-999-6126
		cn.market@abclonal.com.cn
		www.abclonal.com.cn

Product Information

Source Rabbit	Isotype IgG	Purification Affinity purification
-------------------------	-----------------------	--

Storage
Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.