

# Phospho-JunB-S259 Rabbit pAb

Catalog No.: AP0375

1 Publications

## Basic Information

### Observed MW

45kDa

### Calculated MW

36kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human, Mouse, Rat

## Background

Enables sequence-specific double-stranded DNA binding activity. Involved in positive regulation of transcription by RNA polymerase II. Located in nucleoplasm. Part of transcription factor AP-1 complex. Biomarker of Hodgkin's lymphoma and anaplastic large cell lymphoma.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

3726

### Swiss Prot

P17275

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

AP-1; Phospho-JunB-S259

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from HT29 cells, using Phospho-JunB-S259 Rabbit pAb (AP0375).  
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% BSA.