

ABflo® 488 Goat anti-Human IgG (H+L) Ab

Catalog No.: AS181N

Basic Information

Observed MW

Calculated MW

Category

Secondary antibody

Applications

FC

Cross-Reactivity

Human

Conjugate

ABflo® 488. Ex:491nm. Em:516nm.

Background

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Immunogen Information

Gene ID

Swiss Prot

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Goat

Isotype

IgG

Purification

Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.