

Rabbit anti-Mouse MCP-1 mAb (CAP)

Catalog No.: RM17937

Basic Information

Catalog No.
RM17937

Catagory
Elisa Antibody Kit

Application
ELISA, Multiplex

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
/

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Immunogen Information

Immunogen
Recombinant Mouse MCP-1

Cross-Reactivity

No cross-reactivity in ELISA assay with recombinant Human MCP-1


Assay Applications

Mouse MCP-1 Sandwich Immunoassay


		Recommended Concentration	Sample
ELISA	Capture	1-4 µg/mL	Rabbit anti-Mouse MCP-1 mAb (Cat. No. RM17937)
	Detection	0.25-0.0625 µg/mL	Rabbit anti-Mouse MCP-1 mAb (Cat. No. RM17938)
	Standard	7.81-500 pg/mL	Recombinant Mouse CCL2/MCP-1 Protein (Cat. No. RP01607)

Multiplex	Capture	3-20 µg/ml	Rabbit anti-Mouse MCP-1 mAb (Cat. No. RM17975)
	Detection	0.017-2 µg/ml	Rabbit anti-Mouse MCP-1 mAb (Cat. No. RM17937)
	Standard	6.86-5000 pg/mL	Recombinant Mouse CCL2/MCP-1 Protein (Cat. No. RP01607)


Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.



No significant cross-reactivity or interference was observed.



Placed at 37°C for 3 days, the stability of the standard curve all conform to CV <10%.