

Rabbit anti-Human CCL2/MCP-1 mAb

Catalog No.: RM18048

Basic Information

Catalog No.
RM18048

Catagory
Elisa Antibody Kit

Application
Multiplex

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2µm filtered solution in PBS with 0.05% Proclin300, PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

This gene is one of several cytokine genes clustered on the q-arm of chromosome 17. Chemokines are a superfamily of secreted proteins involved in immunoregulatory and inflammatory processes. The superfamily is divided into four subfamilies based on the arrangement of N-terminal cysteine residues of the mature peptide. This chemokine is a member of the CC subfamily which is characterized by two adjacent cysteine residues. This cytokine displays chemotactic activity for monocytes and basophils but not for neutrophils or eosinophils. It has been implicated in the pathogenesis of diseases characterized by monocytic infiltrates, like psoriasis, rheumatoid arthritis and atherosclerosis. It binds to chemokine receptors CCR2 and CCR4. Elevated expression of the encoded protein is associated with severe acute respiratory syndrome coronavirus 2 (SARS‐CoV‐2) infection.

Immunogen Information

Immunogen
Recombinant fusion protein containing a sequence corresponding to amino acids 22-99 of human CCL2/MCP-1 (NP_002973.1).

Cross-Reactivity

Assay Applications

Human CCL2/MCP-1 Sandwich Immunoassay

		Recommended Concentration	Sample
Multiplex	Capture	3-20ug/mL	Rabbit anti-Human CCL2/MCP-1 mAb(Cat. No.RM18047)
	Detection	0.017-2ug/mL	Rabbit anti-Human CCL2/MCP-1 mAb(Cat. No.RM18048)
	Standard	1.37-1000pg/mL	Recombinant Human CCL2/MCP-1 Protein

Validation Data

This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.

