

Rabbit anti-Human NSE/ENO2 mAb

Catalog No.: RMK0225

Basic Information

Catalog No.
RMK0225

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human ENO2 Antibody.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates.

Immunogen Information

Immunogen
Recombinant Human ENO2 Protein

Cross-Reactivity


/

Assay Applications

Human ENO2 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human ENO2 mAb(Cat. No.RMK0225)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human ENO2 mAb(Cat. No.RMK0226)
Standard	6.86-5000pg/mL	Recombinant Human ENO2 Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.