

Rabbit anti-Human TNF RII/TNFRSF1B mAb

Catalog No.: RMK0236

Basic Information

Catalog No.
RMK0236

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human TNF RII/TNFRSF1B Antibody.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

The protein encoded by this gene is a member of the TNF-receptor superfamily. This protein and TNF-receptor 1 form a heterocomplex that mediates the recruitment of two anti-apoptotic proteins, c-IAP1 and c-IAP2, which possess E3 ubiquitin ligase activity. The function of IAPs in TNF-receptor signalling is unknown, however, c-IAP1 is thought to potentiate TNF-induced apoptosis by the ubiquitination and degradation of TNF-receptor-associated factor 2, which mediates anti-apoptotic signals. Knockout studies in mice also suggest a role of this protein in protecting neurons from apoptosis by stimulating antioxidative pathways.

Immunogen Information

Immunogen
Recombinant Human TNF RII/TNFRSF1BProtein

Cross-Reactivity


/

Assay Applications

Human TNFRSF1B multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human TNFRSF1B mAb(Cat. No.RMK0235)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human TNFRSF1B mAb(Cat. No.RMK0236)
Standard	6.86-5000pg/mL	Recombinant Human TNFRSF1B Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.