

Biotinylated Recombinant Human CEACAM8/CD66b Protein

Catalog No.: RP00512B **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	1088	P31997

Tags

C-His&Avi

Synonyms

CD66b; CEACAM8; CGM6; CD67; NCA-95

Product Information

Source	Purification
HEK293 cells	≥ 95 % as determined by SDS-PAGE; ≥ 95 % as determined by HPLC.

Calculated MW	Observed MW
34.3 kDa	60-70 kDa

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from 0.22 μm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

Background

Eosinophils and their products are likely important in the pathophysiology of allergic diseases, such as bronchial asthma, and in host immunity to parasitic organisms. CD66b (CEACAM8, CGM6, NCA-95) is a single chain, GPI-anchored, highly glycosylated protein belonging to the carcinoembryonic Ag supergene family. CD66b is an activation marker for human granulocytes.

Basic Information

Description

Biotinylated Recombinant Human CEACAM8/CD66b Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Gln35-Ser319) of Human CEACAM8/CD66b (Accession #P31997) fused with a C-His&Avi tag at the C-terminus.

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Storage


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

Operational Notes


For your safety and health, please wear a lab coat and disposable gloves for handling.




Validation Data



Biotinylated Recombinant Human CEACAM8/CD66b Protein was determined by Tris-Bis PAGE under reducing conditions.



The purity of Biotinylated Human CEACAM-8 is greater than 95% as determined by SEC-HPLC.



Biotinylated Human CEACAM-8, His Tag immobilized on CM5 Chip can bind Human CEACAM-6, His Tag with an affinity constant of 6.13 μ M as determined in SPR assay (Biacore T200).