

Catalog No.: RP01231 **Recombinant** **3 Publications**

Catalog No.: RP01231

Recombinant

3 Publications

Background

Species	Gene ID	Swiss Prot
Human	9965	O95750

C-His

FGF19

Source	Purification
HEK293 cells	≥ 95 % as determined by SDS-PAGE

Calculated MW	Observed MW
22.09 kDa	25 kDa

< 0.1 EU/μg of the protein by LAL method.

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

 | 400-999-6126

✉ | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Recombinant Human FGF-19 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Phe27-Lys216) of human FGF-19 (Accession #NP_005108.1) fused with a 6×His tag at the C-terminus.


Measured by its binding ability in a functional ELISA. Immobilized Recombinant Human FGF19 at 5 µg/mL (100 µL/well) can bind Recombinant Human FGFR4 with a linear range of 0.9-3.8 µg/mL.

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.
After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.
Avoid repeated freeze/thaw cycles.

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant Human FGF-19 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized Recombinant Human FGF19 at 5µg/mL (100 µL/well) can bind Recombinant Human FGFR4 with a linear range of 0.9-3.8 µg/mL.