

Catalog No.: RP01305 **Recombinant**

Species	Gene ID	Swiss Prot
MERS-CoV	14254594	K9N508

C-mFc

MERS-CoV Spike RBD

Source	Purification
HEK293 cells	≥ 90 % as determined by SDS-PAGE

Calculated MW	Observed MW
51.84 kDa	60-70 kDa

< 0.1 EU/μg of the protein by LAL method.

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

It's been reported that Coronavirus can infect the human respiratory epithelial cells through interaction with the human ACE2 receptor. The spike protein is a large type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which is responsible for recognizing the cell surface receptor. S2 contains basic elements needed for the membrane fusion. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.

Recombinant MERS-CoV Spike RBD Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Glu367-Tyr606) of mers-cov Spike RBD (Accession #YP_009047204.1) fused with a K9N5Q8.


Measured by its binding ability in a functional ELISA. Immobilized Human CD26 (Catalog: RP01300LQ) at 2µg/mL (100µL/well) can bind MERS-CoV Spike RBD with a linear range of 0.5-139.82ng/mL.

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.
After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.
Avoid repeated freeze/thaw cycles.

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant MERS-CoV Spike RBD Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized Human CD26 (Catalog: RP01300LQ) at 2μg/mL (100μL/well) can bind MERS-CoV Spike RBD with a linear range of 0.5-139.82ng/mL.