




Validation Data



Recombinant Human TNFRSF1A/TNF-R1/CD120a Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized recombinant Mouse TNF α at 2 μ g/mL (100 μ L/well) can bind Human TNFRSF1A with a linear range of 4-53 ng/mL.



Recombinant Human TNFRSF1A inhibit TNF- α mediated cytotoxicity in L-929 mouse fibrosarcoma cells in the presence of the metabolic inhibitor actinomycin D. The ED₅₀ for this effect is typically 14-57 ng/mL in the presence of 0.25 ng/mL recombinant human TNF- α .