

Catalog No.: RP01848 **Recombinant**

Species	Gene ID	Swiss Prot
Mouse	14165	O35565

NO-Tag

Fibroblast growth factor 10; FGF-10;
Keratinocyte growth factor 2;Kgf10

Source <i>E. coli</i>	Purification ≥ 95 % as determined by SDS- PAGE.
---------------------------------	---

Calculated MW	Observed MW
19.54 kDa	20-25 kDa

< 1 EU/μg of the protein by LAL method.

Lyophilized from a 0.22 μ m filtered solution of 20mMPB,800mMNaCl pH7.4

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Fibroblast Growth Factor 10 (FGF 10) is an evolutionary conserved secreted growth factor mediating mostly mesenchymal to epithelial signaling. FGF 10 belongs to the FGF 7 subfamily and shares similar biochemical and amino acid sequences with its constituent members (FGF3, FGF 7 and FGF 22). As a paracrine FGF, FGF 10 elicits its biological responses by activating the fibroblast growth factor receptor 2b (FGF R 2b), is crucial for governing proximal distal outgrowth as well as patterning and acts upstream of the known apical ectodermal ridge (AER) marker FGF 8. FGF10 is also implicated in pancreatic cancer, and that overexpression of FGFR2b is associated with metastatic invasion.

Recombinant Mouse FGF-10 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Gln37-Thr209) of Mouse FGF-10 (Accession #NP_032028.1) fused with no tag.

Measured in a cell proliferation assay using 4MBr-5 rhesus monkey epithelial cells. The ED₅₀ for this effect is 3.66-14.64 ng/mL, corresponding to a specific activity of 56.83×10⁴~2.73×10⁵ units/mg.

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.


Avoid repeated freeze/thaw cycles.

For your safety and health, please wear a lab coat and disposable gloves for handling.

Validation Data



Recombinant Mouse FGF-10 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Measured in a cell proliferation assay using 4MBr-5 rhesus monkey epithelial cells. The ED₅₀ for this effect is 3.66-14.64 ng/mL, corresponding to a specific activity of 56.83×10⁴~2.73×10⁵ units/mg.