



Catalog No.: RP02682 **Recombinant**

* □□□□□□□□□□□□□□□□□□□□□□□□


Validation Data



Recombinant Cynomolgus TROP-2/TACSTD2 Protein was determined by Tris-Bis PAGE under reducing conditions.



The purity of Cynomolgus TROP-2 is greater than 95% as determined by SEC-HPLC.



Immobilized Cynomolgus TROP-2 at 0.5µg/ml (100µl/Well). Dose response curve for Anti-TROP-2 Antibody, hFc Tag with the EC₅₀ of 10.2ng/ml determined by ELISA.