

Recombinant Human PGK1 Protein

Catalog No.: RP02987LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	5230	P00558

Tags

C-His

Synonyms

Phosphoglycerate kinase 1;Cell migration-inducing gene 10 protein;Primer recognition protein 2;PGK1;PGKA

Product Information

Source	Purification
HEK293 cells	≥ 95 % as determined by SDS-PAGE.

Calculated MW	Observed MW
45.5 kDa	37-50 kDa

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, pH 8.0.

Reconstitution

Background

Phosphoglycerate kinase 1(PGK1) is an enzyme. It is mainly expressed in spermatogonia and Localized on the principle piece in the sperm. Its expression significantly decreased in the testis of elderly men. PGK1 involved in a critical energy-producing process known as glycolysis. It helps carry out a chemical reaction that converts a molecule called 1,3-diphosphoglycerate, which is produced during the breakdown of glucose, to another molecule called 3-phosphoglycerate during glycolysis. PGK1 may also act as a cofactor for polymerase alpha. The protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions.

Basic Information

Description

Recombinant Human PGK1 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Ser2-Ile417) of human PGK1 (Accession #NP_000282.1) fused with a 6×His tag at the C-terminus.

Bio-Activity

Storage


Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn🌐 | www.abclonal.com.cn

Validation Data



Recombinant Human PGK1 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.