

Catalog No.: RP03186LQ **Recombinant**

Species	Gene ID	Swiss Prot
Mouse	56637	O9WV60

N-His

7330414F15Rik; 8430431H08Rik;
C86142; GSK-3; GSK-3beta; GSK3b

Source	Purification
Baculovirus-Insect Cells	≥ 75 % as determined by SDS-PAGE

49 kDa 45-50 kDa

< 1 EU/μg of the protein by LAL method.

Supplied as sterile 20mM Tris, 500 mM NaCl, 25% glycerol, 0.2 mM DTT, pH 7.4.

Reconstitution

GSK3B is a serine-threonine kinase, belonging to the glycogen synthase kinase subfamily. It contains 1 protein kinase domain, and is expressed in the testis, thymus, prostate, and ovary and weakly expressed in the lung, brain, and kidney. GSK3B is involved in energy metabolism, neuronal cell development, and body pattern formation. Polymorphisms in the GSK3B gene have been implicated in modifying the risk of Parkinson's disease, and studies in mice show that overexpression of this gene may be relevant to the pathogenesis of Alzheimer's disease. GSK3B participates in the Wnt signaling pathway. It is implicated in the hormonal control of several regulatory proteins including glycogen synthase, MYB, and the transcription factor JUN. Phosphorylates JUN at sites proximal to its DNA-binding domain, thereby reducing its affinity for DNA. Phosphorylates MUC1 in breast cancer cells, and decreases the interaction of MUC1 with CTNNB1/beta-catenin. GSK3B also plays an important role in ERBB2-dependent stabilization of microtubules at the cell cortex. It prevents the phosphorylation of APC and CLASP2, allowing its association with the cell membrane. In turn, membrane-bound APC allows the localization of MACF1 to the cell membrane, which is required for microtubule capture and stabilization. GSK3B phosphorylates MACF1 and this phosphorylation inhibits the binding of MACF1 to microtubules which are critical for its role in bulge stem cell migration and skin wound repair. It may be required for early embryo development and neuron differentiation.

Recombinant Mouse GSK-3 beta/GSK3B Protein is produced by Baculovirus-Insect Cells expression system. The target protein is expressed with sequence (Met1-Thr420) of Mouse GSK-3 beta/GSK3B (Accession #Q9WV60) fused with His tag at the N-terminus.

1. The specific activity was determined to be > 20 nmol/min/mg using synthetic Phospho-Glycogen Synthase Peptide-2 (YRRAAVPPSPSLRHSSHPQpSEDEEE) as substrate. |2. Immobilized His-mGSK3B at 10 µg/mL (100 µL/well) can bind biotinylated human HG3C-CTNNB1. EC₅₀ of biotinylated human HG3C-CTNNB1 is 0.15-0.35 µg/mL.


Store at -70°C. This product is stable at $\leq -70^{\circ}\text{C}$ for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Validation Data



Recombinant Mouse GSK-3 beta/GSK3B
Protein was determined by SDS-PAGE under
reducing conditions with Coomassie Blue.