

Recombinant Human Atrophin-1 interacting protein 4/ITCH protein

Catalog No.: RP10137LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human	83737	Q96J02

Tags

NO-Tag

Synonyms

ITCH; E3 ubiquitin-protein ligase Itchy homolog; Itch; EC:2.3.2.26; Atrophin-1-interacting protein 4; AIP4; HECT-type E3 ubiquitin transferase Itchy homolog; NFE2-associated polypeptide 1; NAPP1

Background

ITCH, also called atrophin-1 interacting protein 4 (AIP4), is a HECT-type E3 ubiquitin ligase. It contains an N-terminal Ca

Product Information

Source

E. coli

Purification

≥ 90% as determined by reducing SDS-PAGE.

Calculated MW Observed MW

102.67 KDa

100-140 kDa

Endotoxin

Please contact us for more information.

Formulation

0.22 µm filtered solution of PBS, pH 7.4.

Reconstitution

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn