

Recombinant Human Ubiquitin(K0) protein

Catalog No.: RP10176LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Human		P0CG47

Tags

No Tag

Synonyms

Ubiquitin

Product Information

Source	Purification
<i>E. coli</i>	≥ 95 % as determined by SDS-PAGE.

Calculated MW **Observed MW**

Endotoxin

Please contact us for more information.

Formulation

20 mM Tris, 150 mM NaCl, 2 mM βME,
10% Glycerol

Reconstitution

Background

All lysine residues in Ub(K0) were substituted with arginine residues, thereby Ub(K0) can only support monoubiquitination, but not polyubiquitination.

Basic Information

Description

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.


Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Validation Data



Recombinant Human Ubiquitin(K0) protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.